Varolast® Plus

Zinkleimbinde, längselastisch

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EN Instructions for use

General description of the device

Varolast Plus are ready-to-use lengthwise elastic zinc paste bandages. The bandages have a carriage material made of viscose and polyamide, and use zinc paste as impregnation mass, which contains zinc oxide, glycerin and gelatin. Zinc oxide is added due to technical reasons to achieve a semi-rigid consistency of the bandages; an ancillary pharmacological action is not intended. The bandages are white in colour. The non-active devices are available in various widths and lengths.

Intended purpose and indications

Varolast Plus are non-sterile, single-use zinc-paste bandages used by healthcare professionals on intact skin of humans (above 3 years of age) in clinical and home environment. They can be used for long-term compression therapy (ABPI 0.9 – 1.3) (decongestion phase) of:

- varicose veins (e.g. varicose in pregnancy, supportive in the invasive treatment of varicosis)
- chronic venous insufficiency (CVI) with or without venous leg ulcers (VLU) (e.g. CVI according to CEAP classification C3 – C6, primary and secondary prevention of leg ulcer venosum, venous insufficiency and angiodysplasia)
- thrombosis/thrombophlebitis (e.g. thrombophlebitis as well as condition after healed thrombophlebitis, deep leg and arm vein thrombosis, condition after deep vein thrombosis, post-thrombotic syndrome, thromboprophylaxis)
- chronic oedemata of other origin (e.g. oedema in pregnancy, post-traumatic oedema, hormonal oedema, lipodema, stasis conditions due to immobility (arthrogenic congestion syndrome, paresis), medically induced oedema (e.g. calcium antagonists, isosorbide dinitrate, lithium ointment, sex hormones))
- lymph oedema

Further, they can be used as a semi-rigid support bandages. The compression bandages can also be used in combination with primary and secondary dressings.

Composition

Bandage: 71% viscose, 29% polyamide; Zinc paste: Water, zinc oxide, glycerin, and gelatin - packed in an aluminum coated foil pouch.

Contraindications

Varolast Plus must not be used in case of:

- advanced peripheral arterial occlusive disease (PAOD)
- decompensated heart failure (NYHA III + IV)
- phlegmasia coerulea dolens
- ankle-brachial pressure index (ABPI) > 1.3 and ≤ 0.5
- septic phlebitis
- florid erysipelas

Precautions

Special precautionary measures are required in the event of:

- pronounced weeping dermatoses
- pronounced sensibility disorders of the extremities (e.g. in the case of advanced diabetic peripheral polyneuropathy)
- primary chronic rheumatoid arthritis
- incompatibility with compression bandage material or zinc paste ingredients

The application technique, e.g. with regards to stretch and/or padding, may have to be adapted to the patient's anatomical and pathophysiological characteristics, especially in case of fragile or aged skin.

If the foot pulse is not palpable, a Doppler ultrasound should be carried out to determine the ABPI prior to starting therapy. The right measurement method must be considered for the ABPI. During the use of the compression bandages, the ABPI must be controlled regularly. The therapy must be adapted accordingly.

If an inelastic system is created by close to full stretch application, high ankle pressures should be closely monitored by a physician.

The devices should not be used on children below the age of 3 years.

When applying Varolast Plus on children (from 3 years on), the pressure should be reduced appropriately. If at any time during therapy the patient develops severe pain, numbness, tingling or noticeable discolouration of the toes, a healthcare professional must be contacted immediately, and the compression bandages must be removed.

Further remarks

- Do not let the bandages get wet during showering/bathing.
- Do not use the bandages as a primary dressing or directly on open wounds.
- Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request.

Incident reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation (EU) 2017/745 on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Product disposal

In order to minimize the risk of potential infection hazards, or environmental pollution, disposable components of the medical device should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards. Dispose of the medical device with regular hospital waste.

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Paul Hartmann AG Paul-Hartmann-Straße 12 89522 HEIDENHEIM, GERMANY www.hartmann.info