

ReadyWrap Foot CT

Compression garment

REF 141726 - 141765

May 2022

1. Composition of the product

The ReadyWrap Foot CT compression garment consist of:

- Nylon
- Polyurethane
- Elastane

This product data sheet covers the following item:

REF	Product designation	Version	Size	Color
141726	ReadyWrap Foot CT	left	REG S	beige
141727	ReadyWrap Foot CT	right	REG S	beige
141728	ReadyWrap Foot CT	left	REG S	black
141729	ReadyWrap Foot CT	right	REG S	black
141746	ReadyWrap Foot CT	left	LONG S	beige
141747	ReadyWrap Foot CT	right	LONG S	beige
141748	ReadyWrap Foot CT	left	LONG S	black
141749	ReadyWrap Foot CT	right	LONG S	black
141730	ReadyWrap Foot CT	left	REG M	beige
141731	ReadyWrap Foot CT	right	REG M	beige
141732	ReadyWrap Foot CT	left	REG M	black
141733	ReadyWrap Foot CT	right	REG M	black
141750	ReadyWrap Foot CT	left	LONG M	beige
141751	ReadyWrap Foot CT	right	LONG M	beige
141752	ReadyWrap Foot CT	left	LONG M	black
141753	ReadyWrap Foot CT	right	LONG M	black
141734	ReadyWrap Foot CT	left	REG L	beige
141735	ReadyWrap Foot CT	right	REG L	beige
141736	ReadyWrap Foot CT	left	REG L	black
141737	ReadyWrap Foot CT	right	REG L	black
141754	ReadyWrap Foot CT	left	LONG L	beige
141755	ReadyWrap Foot CT	right	LONG L	beige
141756	ReadyWrap Foot CT	left	LONG L	black
141757	ReadyWrap Foot CT	right	LONG L	black
141738	ReadyWrap Foot CT	left	REG XL	beige
141739	ReadyWrap Foot CT	right	REG XL	beige
141740	ReadyWrap Foot CT	left	REG XL	black
141741	ReadyWrap Foot CT	right	REG XL	black

REF	Product designation	Version	Size	Color
141758	ReadyWrap Foot CT	left	LONG XL	beige
141759	ReadyWrap Foot CT	right	LONG XL	beige
141760	ReadyWrap Foot CT	left	LONG XL	black
141761	ReadyWrap Foot CT	right	LONG XL	black
141742	ReadyWrap Foot CT	left	REG XXL	beige
141743	ReadyWrap Foot CT	right	REG XXL	beige
141744	ReadyWrap Foot CT	left	REG XXL	black
141745	ReadyWrap Foot CT	right	REG XXL	black
141762	ReadyWrap Foot CT	left	LONG XXL	beige
141763	ReadyWrap Foot CT	right	LONG XXL	beige
141764	ReadyWrap Foot CT	left	LONG XXL	black
141765	ReadyWrap Foot CT	right	LONG XXL	black

2. Packaging, structure and composition

2.1 Unit container

- See. 2.2

2.2 Shelf Container

- 1 piece = 1 PE flat bag with pressure cap
- 1 instructions for use consisting of cellulose
- 1 liner

2.3 Transit Container

- 10 pieces = 1 folding boxes
consisting of corrugated cardboard box (cellulose)

3. Manufacture

ReadyWrap Foot CT are produced according to specification in hygienic conditions and packed as described in its relevant packaging specification.

4. Description

ReadyWrap Foot CT adjustable compression garments with Velcro brand fasteners are designed as an easy to use compression therapy for patients with venous and lymphatic disease.

5. Properties

- easy to use and apply with three-step colour-coded hook-and-loop fasteners
- short stretch properties - therapeutic working pressure and well tolerated lower resting pressure

- ReadyWrap Foot CT edges overlap to ensure that the legs or arms are being completely covered
- available in black or beige

6. Intended purpose

Compression garment for phlebological, lymphological indications i.e. venous disease, chronic oedema and lymphoedema. To maintain or reduce fluid volume and shape of upper and lower extremity limbs.

7. Medical device classification

ReadyWrap Foot CT are medical devices of class I according to rule 1. (according to regulation (EU) 2017/745)

8. Biological evaluation and biocompatibility

The materials in general in ReadyWrap Foot CT are used in the compression area since years and are safe if the product is used on intact skin appropriately.

There is no known risk involved in the use of the medical device ReadyWrap Foot CT and it is designed, manufactured and packaged in such a way that it will not compromise the clinical condition or the safety of patients, or the safety and health of users and other persons when used under the conditions and for the purposes intended.

9. Requirements regarding the safety/harmlessness of the constituents

Currently we do not have any information on questionable ingredients.

10. Stability

Stored appropriately (store in a dry place, protected from dust and light), ReadyWrap Foot CT have a shelf life of 3 years (ReadyWrap garment: 5 years, Liner: 3 years).

11. Disposal

The user is advised to observe current national legislation, norms and guidelines, regulating the disposal of medical refuse.

Packaging materials must also be disposed of in compliance with applicable national requirements.

Lohmann & Rauscher International GmbH & Co. KG
D-56579 Rengsdorf
signed by

Peter Grillitsch
Director Global Portfolio Management